



TECNO TREND

Aseptic Powder Machines

MICRO POWDER DOSING MACHINES (WITHOUT ROTATING DISK)			
Technical data	FY100T	FY200T	FY300T
Output – single dose (vials/hr)	6.000	12.000	18.000
Vial dimensions diam.	8 mm - 60 mm	8 mm – 28 mm	8 mm - 23 mm
Vial dimensions h	40 mm -160 mm	40 mm -160 mm	40 mm - 160 mm
Minimum dose (Min. dosages can change according to the powder specific weight)	0,02 cc	0,02 cc	0,02 cc
Maximun dose (Max. dosages can change according to the powder specific weight)	7,5 cc (1,5 g. with single dosage)	7,5 cc (1,5 g. with single dosage)	7,5 cc (1,5 g. with single dosage)
Power consumption	8 Kw	8 Kw	8 Kw
Compressed air consumption	100 NI/min	100 NI/min	100 NI/min
Compressed air pressure in bar	6	6	6
Vacuum pump	100m3/h	100m3/h	100m3/h
PLC standard	SIEMENS	SIEMENS	SIEMENS
Net weight (kg)	1100	1100	1200
Gross weight (kg)	1400	1400	1500